



Ion Detoxification

Operating Instructions

Please read the manual in detail before using

Operating Instructions

A. Product Principle:

The water is decomposed by this ion cleanse device into large numbers of energy ions (positive ions and negative ions), which can strengthen the vigor of body cells, promote the metabolism and improve the sleep quality so as to furbish the health gradually.

B. Patients and practitioners of the system have reported many great effects from using the system. They include:

1. Increased joint movement
2. Joint pain relieve from arthritis
3. Better sleep
4. A better well being
5. Increased energy

C. Operation

1. To begin, plug in the power cord, the array, and the wrist strap into the machine. Next, put the array into the wash pan filled with water, and put on the wrist strap. Finally, turn on the power.
2. When you turn on the power switch, the screen shown is as follows

```
Enter Time 10-30  
Time=30 V3.0
```

Default time is Time=30 minutes. To adjust for more or less time, use the (+) or (-) key.

```
Time=06:28 pos+  
T=30 M=3 E=1.0
```

Time=06:28 pos+ positive polarity, pos- negative polarity
T=30 minutes time, M=3 mode three, E=1.0 current output, The instrument automatic turns off when the timer ends.

3. During initial operation, mode 3 is default. The current output is selected by choosing a different mode. The current output is as follows: Mode 1=0.6 amps, Mode 2=0.8 amps, Mode 3=1.0 amps. To select a different Mode, press (-) or (+) key. Pressing (ESC) will return to menu. Because there are differences in local

water quality, it may be necessary to add a pinch of salt to the water to maintain the desired current output.

4. This product has the recorded time and the number of times function. When you need to examine the record, simultaneously press (1) and (ENT) on the keypad, and turn on the power switch. this time the screen reads:

```
00298 Mins
93 Time
```

The first line shows total operating time, the second line shows the total numbers of times the machine has been used. When the machine has achieved 10 minutes time, it will counts as one use. Use (ESC) to go back to the menu.

5. To delete the time and counter, press (2) and (ENT). and turn on the power switch. This time the screen reveals:

```
PAS: V3.0
Clear Contents
```

The default password is “666666” then press (ENT) The display reads:

```
Process: OK
Mins=00000
```

Press (ESC) to return to the menu.

6. To make your own password, press (3) and (ENT). and turn on the power switch. The display reads:

```
PAS: V3.0
Change Password
```

First, input the original six numbers, (666666) and press (ENT), this time the screen reveals:

```
Process: OK
New PAS: 123456
```

Now input your new password, a six number password, and press (ESC). Remember your new password. The display should be the menu.

7. In the following cases, the machine will make sounds or an error character will show on the screen.
- The array is not in the water or the array is broken.
 - The wrist strap is not wrapped in the proper way or the wrist strap is broken.
 - The current power is more than 2.5 amps.
8. At any time, you may use the bamboo charcoal, far infrared waistband. Connect the waistband to the machine, put the waistband tightly around your belly, turn on the infrared power switch, choose mode L(low), M(medium), or H(high) to start far infrared massage.

D.Attention

1. Do not touch power while your hand is wet, please pay attention to your personal safety.
2. Avoid striking the machine.
3. Use clean cloth to polish the dust on the machine's surface, do not use any chemicals.
4. Be gentle when pressing the keys to protect the machine under normal use.
5. You must prop up the machine with the handle while the machine is operating.

E.Packing list

Ion Detoxification	Array	Power Cord	Wrist strap	Aluminum Storage Box	Far Infrared Waist Belt	Instructions
1	1	1	1	1	1	1

F.Operating instructions

1. Use ordinary tap water. Do not use distilled water. All water quality is different. Only add a little salt at a time. Adding salt will promote the Ion Detoxification System to conduct electricity. The capability of conducting mainly lies in the waters conductivity quality. The capability of the conducting of the water determines the ion effects and impacts the function of detoxification. If the ampere readout is 1.5 amperes or below, add a little salt to enhance the capability of the conducting. But, for maintaining the machine life, please avoid adding salt.
2. Adults need 30 minutes per section.
3. 8-15 year old children must use with adult supervision for no more than 15 minutes.
4. Under 50 years, OK to use it every other day for 14 continuous treatments, then stop for three weeks before starting the 14 treatment cycle. Over 50 years, OK to use it every three days for 14 continuous treatments, and then stop for four weeks before starting the 14 treatment cycle. Do not use more than three days in a row. Do not use more than twice per day.
5. The array is needs to be cleaned after every use by soaking in a combination of 30 parts water and 1 part of phosphoric acid or household bleach. Soak at least 5 minutes.

G.Precautions

- 1.People who have battery operated implants(e.g..pace makers)should never use
- 2.Anyone who has received an organ transplant should never use.
- 3.Pregnant women or nursing woman should never use.
- 4.The exterior wounds,if having,can not dip in the water when ion cleansing.
- 5.Less than 8 years old child should not use.
- 6.People who have hypoglycemia should have food before use.

The Far Infrared Ray used with The Bamboo Charcoal Waist Belt Instructions

The latest studies have shown that bamboo charcoal made from 5 year old bamboo from remote mountains, refined through high temperature have a great absorption effect to purify the air or water and adjust humidity. The bamboo charcoal also emits far infrared ray, negative ions, and absorbs electromagnetic waves. The waist belt made of bamboo charcoal can warm your body, relieve muscle aching, increase the cells' activities, and accelerate the blood circulation.

Far Infrared Ray used with the bamboo charcoal waist belt can help with

- 1.Decreases joint stiffness.
- 2.Relieves muscle spasms.
- 3.Increases blood flow.
- 4.Leads to pain relief.
- 5.Affects soft tissue injury.
- 6.Increases the extensibility of collagen tissue.
- 7.Assists in resolution of inflammatory infiltrated, edema, and exudes.

The most light we see every day is actually a combination of many different types. All lights are electro-magnetic waves. The wavelengths of lights determine their properties. Lights could be visible, or invisible. Visible lights could also present different colors. As shown in above image, the Sun lights spectrum ranges from ultraviolet and violet on the short wavelength side, to red and infrared on the long wavelength side.

Different types of light interact with our body differently. In general, the shorter the wavelength of a light, the harmful it is to living cells. Lights, like ultraviolet, are very harmful to living cells. The most beneficial lights are those whose wave

lengths are near the red light. Far Infrared Therapy, or FIR therapy or Infrared Therapy uses far infrared light which is part of the Sun light spectrum whose wave lengths fall just below (what the word infra means) the visible red light.

Although the wavelengths of FIR are too long for the eyes to perceive, we can experience its energy as gentle, radiant heat, which can penetrate up to [3.5] inches beneath the skin. So Far Infrared Therapy, or FIR Therapy is also a kind of light therapy and heat therapy, it works to enhance health through following principles:

- Improve micro circulation by exerting strong rotational and vibration effects in the molecular level, which are biologically beneficial.

- Increase the delivery of oxygen and nutrients in the blood cell to the body's soft tissue areas.

- Enhance white blood cell function, thereby increasing immune response and the elimination of foreign pathogens and cellular waste products.

- Remove accumulated toxins by improving lymph circulation which are often at the core of many health problems. The build up of toxins in a healthy body could block the normal blood circulation and impair the cellular energy. When FIR waves are applied, the water molecules that encapsulates the toxins get heat up, and start to vibrate. This vibration reduces the ion bonds of the atoms that are holding together the molecules of water. As the breakdown of the water molecules occurs, encapsulated gases and other toxic materials are released and the body gets rejuvenated.

- Promote regeneration and fast healing.

- Increase metabolism between blood and tissue.

Stimulate the hypothalamus, which controls the production of neurochemicals involved in such biological processes as sleep, mood, pain sensations, and blood pressure

Through above mechanism, Far Infrared Therapy with FIR therapeutic clothing products could put the body under far infrared treatment any time of the day, and however long the user wants. With low dose of FIR energy, there is no risk of burning of the body cell. Far Infrared Therapy has been shown to have good effects in treating following symptoms:

- Arthritis
- Muscle pain and spasm
- Shoulder or joint stiffness
- Raynaud's syndrome

Far Infrared Therapy has also been shown to aid the process of weight loss, promote the growth of new hair, improve immune system and the overall health.

Far infrared therapy follows the holistic healing principle. It does not only treat a symptom, it enhances the whole body function, the treatment of the symptom simply comes along naturally.

The system does not claim to cure or treat any disease or injury. It assists your body to re-balance its bio-energy fields, and stimulates the body for self-detoxification. The body's organs will naturally function better when the electro-magnetic fields are charged. The ion energy charged footbath will help you feel young again.

Warranty Information Card

This warranty covers defects in material and workmanship for one year from the original date of delivery. We will replace any part which in our opinion is defective; if the housing has not been subjected to tampering, alteration, or improper use after delivery and has not been repaired by anyone except us. Our obligation does not include the cost of transportation. We are not responsible for damage in transit, and claims for such damage should be presented to the carrier by the customer. The housing is not warranted against misuse or use in abnormal temperature conditions.

For servicing under this warranty, return any defective part to YOUR DISTRIBUTOR within the period listed below.

Warranty Registration Card

To register your warranty for this product, complete the information below.

Name

Address

City/State/Zip

Telephone

Date Purchased



Color or Particle	Material or Area of the Body Being Detoxified
Black	Detoxifying from liver
Black Flecks	Heavy metals
Brown	Detoxifying from liver, cellular debris, tobacco
Dark Green	Detoxifying from gallbladder
Orange	Detoxifying from joints
Red Flecks	Blood clot material
White Cheese-Like Particles	Most likely yeast
White Foam	Detoxifying from Lymphatic system
Yellow-Green	Detoxifying from kidneys, bladder, urinary tract, female/prostate area

MIRACLE CO.,LTD

Address: Room B601 Hengsheng Building B , No. 205 Huanshi Zhong Road Guangzhou
China /POST NO.:510010

TEL:86-20-86673107 FAX:86-20-86672750

Email: admin@miracle-china.com

Skype: miracle-china

MSN:china.miracle@hotmail.com QQ:286372662

<http://www.miracle-china.cn>

www.miracle-china.cn