

Function Display of the Control Panel




Names of parts



Main energy cushion (extendable)

Every function of the rehabilitation device is exported through this energy cushion. The user sits or lies on the energy cushion. Via the extension jack, $2\sim4$ more energy cushions can be linked together.



Auxiliary energy cushion (non-extendable)

Every function of the rehabilitation device is exported through this energy cushion. The user sits or lies on the energy cushion, which is only for one person.




Energy cushion specially for sleep field (non-standard configuration, separately purchase is needed)

It is for sleep field exportation. Receiving treatments during sleep is both time-saving and more efficient.



Microprobe (pointed type)

The microprobe (an electronic pen), via high-voltage discharge, produces air ions flow, high-voltage micro current and local stimulation effect, thus can lower the excitability of some nerves, promote blood circulation and remove blood stasis, diminish inflammation and ease pain and quickly eliminate symptoms of local numbness, swelling and pains.



Isolation cushion

While in treatment, it is placed on the surface, separating the user and the earth to ensure the safety of treatment.

Technical Specifications

Field Energy Mode	High Potential/ Negative Potential	Sleep Field	High-cycle
Field Energy Strength	$\pm 1500 m V \ \pm 3500 m V \ \pm 5000 m V \ + 29000 m V$	1500V	2000 V
Operating Time Setting	30 minutes ~6 hours		
Dimensions	Ф25cm* h50cm		
Working Environment Temperature	5~40℃		
Relative Humidity	≤80%		
Power Consumption	60W~120W		
Rated Power Supply	AC220V±220V,50H	Hz±1Hz	

MIRACLE CO.,LTD

Address: Room B601 Hengsheng Building B, No. 205 Huanshi Zhong Road Guangzhou China /POST NO.:510010 TEL:86-20-86673107 FAX:86-20-86672750

Email: admin@miracle-china.com

Skype: miracle-china

atmail.con.
inatch MŚN:china.miracle@hotmail.com QQ:286372662

http://www.miracle-china.cn